

# 700\*230mm

EN  
fluorecare®

CE 2934

## PRODUCT NAME

### SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay)

REF MF-71

## WHAT DOES THE KIT TEST?

The fluorescent SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit is applicable to the simultaneous qualitative detection and differentiation in vitro of novel Coronavirus (SARS-CoV-2) antigen, Influenza A virus antigen, Influenza B virus Antigen, and RSV antigen in nasal swab samples.

If it is used for COVID-19 test, it can detect COVID-19 and the infectious disease COVID-19 caused by SARS-CoV-2 in symptomatic patients within 7 days of onset and in the infectious disease COVID-19 caused by Influenza A/B or RSV.

For a direct diagnosis use only. For self-testing:

User age requirement

This kit is suitable for users aged 14 years old and above.

People aged 14-24 should not use the kit by themselves. Sample collection and testing with this kit should be performed by a medical professional.

People aged 14-17 may use this kit to collect samples and test samples under the supervision of adults or parents aged 14-17 years old.

The person supervising use of the kit should be at least 18 years old and have received training on how to use the kit.

For over 18 years old, if it is recommended that family members or caregivers (aged 18-60 years old) perform sample collection and testing of nasal swab samples.

## BACKGROUND

The novel coronavirus belongs to the genus *Beta-CoV*. It is an acutely respiratory infectious disease. People are generally asymptomatic, but may also be asymptomatic. Based on the current epidemiological investigation, the incubation period is 1-14 days, and most 3-7 days. The main manifestations include fever, dry cough, fatigue, and sore throat.

Influenza is a contagious respiratory illness caused by influenza viruses. Some people, such as the elderly, children, and those with chronic diseases, are at higher risk of developing severe illness.

Severe illness may result in hospitalization or death. Some specific populations, such as the elderly, children, and those with chronic diseases, are at higher risk of developing severe illness.

Influenza viruses regularly spread in people, and are responsible for seasonal flu each year. Influenza viruses can be spread through droplets from the nose and mouth.

Respiratory syncytial virus (RSV) belongs to the Paramyxoviridae family of viruses.

It is a common respiratory tract infection in young children, especially in the winter months.

RSV can cause serious complications, such as bronchitis, pneumonia, in infants under 6 months of age, and upper respiratory tract infections, such as rhinitis and the common cold, in older children and adults, and bronchitis in the elderly.

## PRINCIPLE

The SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit is a colloidal gold chromatographic immunoassay. The SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit in the collected sample will be combined with the SARS-CoV-2 (or Influenza A/B or RSV) antibody-coated colloidal gold particles. Due to chromatography, the SARS-CoV-2 (or Influenza A/B or RSV) antibody-coated colloidal gold complex along the microtunnel membrane. Within the detection line, there is a color reaction, showing a purple band. When colloidal gold-labeled SARS-CoV-2 or Influenza A/B or RSV antigen binds to the corresponding antibody-coated colloidal gold, a color reaction occurs.

Anti-mouse IgG is red or brown. When the reaction is over, the results can be interpreted by visual observation.

## SPECIFICATIONS

1 box, 2 test boxes, 5 tests.

Components

[REF] MF-71a [REF] MF-71b [REF] MF-71c

Test Card (including the desiccant)

1 cassette 2 cassettes 3 cassettes

Constituent: Plastic shell, desiccant, the test strip, Type A (Influenza A/B antibody), Type B (Influenza A/B antibody), and the conjugate pad containing colloidal gold particles. The test strip contains Type A (Influenza A/B antibody) and Type B (Influenza A/B antibody). Other components include PVC pad and absorbent paper.

Instructions for Use 1 copy 1 copy 1 copy

Sterile nasal swabs 1 piece 2 pieces 5 pieces

Profiled Sample treatment tube 1 tube 2 tubes 5 tubes

Sample treatment solution (STS) 1ml sterile saline solution 0.5ml / tube

MATERIALS REQUIRED BUT NOT PROVIDED

Time

## STORAGE CONDITIONS AND EXPIRATION DATE

1. Store the kit at 2-30°C, in a dry place, protected from light. Test is valid for 24 months.

2. The Test Card must remain in the sealed pouch until use. Once the test card is opened, the test should be performed within 2 hours.

## TEST PROCEDURE

Use a disinfected to disinfect your hands before testing.

Close your nostrils so that the test will be performed.

Before testing, read the instructions for use carefully, and bring the testing kit and samples to room temperature (20-25°C). The test should be performed at 20-25°C. If the kit is removed from the refrigerator, allow it to stand

at room temperature (20-25°C) for 4 minutes before testing.

1. Wash your hands with the soap provided and use the blue stopper. The purpose of the blue stopper is to prevent liquid from leaking during transportation. The blue stopper should be removed before use. Insert the treatment tube into the hole of the kit or use other items to hold the treatment tube upright.

2. Turn open the fast tab, take out the test card, and use it as soon as possible within 1 hour.

3. Nasal swab collection method:

1) Carefully pull the nasal swab from the packaging. Avoid touching the end of the nasal swab with the nose.

2) Rotate the nasal swab on the moist wall (moisture number) 5 times to ensure adequate sampling.

3) Repeat the process in the right nostril with the same nasal swab, collecting from both nostrils to ensure an adequate sample.

4. Place the swab sample into the tube, then break the swab at the break line and leave the lower half in the treatment tube. Secure the cap back tightly.

5. Squeeze the swab 10 times in the test tube for 1 minute of sample reaction time. Unscrew the tip cover of the tube, tip the cover on the swab and do not touch the swab inside the tube. After drying the sample, the tip of the swab is dry and the swab inside the tube is wet.

6. The test is kept at room temperature for 12 minutes to observe the test results. Results observed after 12 minutes are invalid. While waiting for results, DO NOT touch the test card or move the test card from its original position.

7. The test is valid for 24 hours.

8. Negative result: No color change in the control line (C) and test line (T).

9. Positive result: A color change in the test line (T) is observed.

10. Invalid result: The test card is damaged or the test card is not valid.

11. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

12. The result should be read at 12 minutes. Please do not read the result after 20 minutes.

13. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

14. This kit is for a single use only. Do not reuse any component of the test kit.

15. Negative result does not mean that you are not infected. You may also be infected but not detected.

16. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

17. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

18. The test is valid for 24 hours.

19. Negative result: No color change in the control line (C) and test line (T).

20. Positive result: A color change in the test line (T) is observed.

21. Invalid result: The test card is damaged or the test card is not valid.

22. The test is valid for 24 hours.

23. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

24. This kit is for a single use only. Do not reuse any component of the test kit.

25. Negative result does not mean that you are not infected. You may also be infected but not detected.

26. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

27. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

28. The test is valid for 24 hours.

29. Negative result: No color change in the control line (C) and test line (T).

30. Positive result: A color change in the test line (T) is observed.

31. Invalid result: The test card is damaged or the test card is not valid.

32. The test is valid for 24 hours.

33. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

34. This kit is for a single use only. Do not reuse any component of the test kit.

35. Negative result does not mean that you are not infected. You may also be infected but not detected.

36. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

37. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

38. The test is valid for 24 hours.

39. Negative result: No color change in the control line (C) and test line (T).

40. Positive result: A color change in the test line (T) is observed.

41. Invalid result: The test card is damaged or the test card is not valid.

42. The test is valid for 24 hours.

43. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

44. This kit is for a single use only. Do not reuse any component of the test kit.

45. Negative result does not mean that you are not infected. You may also be infected but not detected.

46. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

47. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

48. The test is valid for 24 hours.

49. Negative result: No color change in the control line (C) and test line (T).

50. Positive result: A color change in the test line (T) is observed.

51. Invalid result: The test card is damaged or the test card is not valid.

52. The test is valid for 24 hours.

53. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

54. This kit is for a single use only. Do not reuse any component of the test kit.

55. Negative result does not mean that you are not infected. You may also be infected but not detected.

56. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

57. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

58. The test is valid for 24 hours.

59. Negative result: No color change in the control line (C) and test line (T).

60. Positive result: A color change in the test line (T) is observed.

61. Invalid result: The test card is damaged or the test card is not valid.

62. The test is valid for 24 hours.

63. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

64. This kit is for a single use only. Do not reuse any component of the test kit.

65. Negative result does not mean that you are not infected. You may also be infected but not detected.

66. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

67. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

68. The test is valid for 24 hours.

69. Negative result: No color change in the control line (C) and test line (T).

70. Positive result: A color change in the test line (T) is observed.

71. Invalid result: The test card is damaged or the test card is not valid.

72. The test is valid for 24 hours.

73. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

74. This kit is for a single use only. Do not reuse any component of the test kit.

75. Negative result does not mean that you are not infected. You may also be infected but not detected.

76. Positive result does not mean that you are definitely infected. You may also be infected but not detected.

77. If the test is invalid, repeat the test procedure and the test should be repeated. If the result is still invalid, please contact customer service.

78. The test is valid for 24 hours.

79. Negative result: No color change in the control line (C) and test line (T).

80. Positive result: A color change in the test line (T) is observed.

81. Invalid result: The test card is damaged or the test card is not valid.

82. The test is valid for 24 hours.

83. Various components of different batches of reagents cannot be used interchangeably in order to avoid inaccurate results.

84. This kit is for a single use only. Do not reuse any component of the test kit.

85. Negative result does not mean that you are not infected. You may also be infected but not detected.